

# Model M Quick Start Guide



2 Insert freshly charged battery (DLI88)



**IMPORTANT!**


Slide battery in from the side. DO NOT press in from top to avoid breaking the gold spring pin contacts:




6 When GREEN LED is flashing, unit is operating correctly



4 Wait 2 seconds for screen to clear  
It will sweep all white then black



3 Set switch to "battery" position (down)



5 Press RED button 3 times firmly

Visit [micro-pulse.com](http://micro-pulse.com), go to Products, and click on the model M1 to find links to documents and instructional videos

# NOT FDA APPROVED

## Intended for experimental purposes only

### Control Panel Overview:

ICES model M series  
 Rev: M-20170706  
 (C) 2017 Micro-Pulse  
 < SET >

POWER LEVEL = 9  
 [ ■■■■■■■■ - - - - - ] 9  
 -- B5 – C5 --  
 < CHANGE > < SET >

SELECT PROTOCOL  
 -- B5 – C5 --  
 < CHANGE > < SET >

-- B5 – C5 --  
 [ ■■■■■■■■ - - - - - ] 9  
 - to reset: power OFF -

**Powering ON display**  
 Press RED button to start your previous program immediately

**Adjust Power Level display**  
 Press <WHITE> to adjust power: 15 (high) to 1 (low), then back to 15  
 Press <RED> to begin immediately

**Select Protocol display**  
 Press <WHITE> to select protocol: See all protocols in list at right →→  
 Press <RED> to begin immediately

**System Running display**  
 Displays Protocol and Power Level  
 Must switch power OFF then ON to change power or protocol settings

Green flashing = pulses being sent  
 Red steady = ERROR or malfunction  
 Yellow = REST (no pulses during rest)



**LEDs:**

For general use:  
 Experiment on pain and inflammation

- B5 – C5 --
- A9 --
- P2 -- (SomaPulse, AllevaWave, ...)
- Omni 8 --

Experiment with Schumann frequencies

- Schumann 1 -- (7.83 pps)
- Schumann 2 -- (7.83, 14.3 pps)
- Schumann 3 -- (7.83, 14.3, 20.8 pps)
- Schumann 4 -- (7.83, 14.3, 20.8, 27.3 pps)
- Schumann 5 -- (7.83, 14.3, 20.8, 27.3, 33.8 pps)

Experiment with fixed constant frequencies

- 1 pps -- continuous bipolar pulses at 1 pulse per second
- 2 pps -- continuous bipolar pulses at 2 pulses per second
- 3 pps -- continuous bipolar pulses at 3 pulses per second
- 4 pps -- continuous bipolar pulses at 4 pulses per second
- 5 pps -- continuous bipolar pulses at 5 pulses per second
- 10 pps -- continuous bipolar pulses at 10 pulses per second

Experiment with low-power TMS

- scTMS 10pps 30 minutes -- Low power TMS simulator
- scTMS 10pps 60 minutes -- Low power TMS simulator

Experiment with brainwave entrainment

- alpha wave -- 10 Hz to 13 Hz (10 minute cycle time)
- beta1 wave -- LOW: 12.5 Hz to 16 Hz (4 minute cycle)
- beta2 wave -- MED: 16.5 Hz to 20 Hz (4 minute cycle)
- beta3 wave -- HIGH: 20 Hz to 28 Hz (4 minute cycle)
- delta wave -- 1.5 Hz to 3 Hz (20 minute cycle time)
- theta wave -- 5 Hz to 6 Hz (20 minute cycle time)
- mu wave -- 8.5 Hz to 11 Hz (10 minute cycle time)
- SMA wave -- 13 Hz to 15 Hz (10 minute cycle time)
- gamma wave -- 32 Hz to 47 Hz (10 minute cycle time)

Experiment using rest periods between cycles

- B5 – C5 -- REST 5 -- 5 minute rest at end of each cycle
- A9 -- REST 5 -- 5 minute rest at end of each cycle
- P2 -- REST 5 -- 5 minute rest at end of each cycle
- Omni 8 -- REST 5 -- 5 minute rest at end of each cycle